



CUSTOMER					SHIP TO			
Company					Name			
Account #			Date		Address			
O Quote	o P0#				City		State	
Buyer					Post Code	Phone		
Mark For					Email			
O Quote has Additional Sunrise Products with different Quote / Order #			Quote / Order #					

Technical Data Overview			
Maximum User Weight	120 Kg (265 lbs)	Motor	150 W
Wheelchair Rear Wheel Sizes	24" (540)	Operating Voltage	24V
Total Device Weight	12 kg (26.4 lbs)	Maximum Motor Noise	60Db
Per Battery Weight (4.0Ah)	0.6 kg (1.3 lbs)	Battery Type	Lithium Ion
Optimal Range Per Charge	15 km (9.3 mi)	Device Charger Type	28V DC 5.7A 160W
Maximum Speed	6 kph (3.7 mph)	Docking Charger	28V DC 1.5A IP44
Maximum Slope	6°	Control Box	3.7V 1Ah 3.7Wh
Motor Wheel Size	24" (540) Tires	Bluetooth Standard	BLE (125Kbps Coded PHY)

BASE			
a OPTION	Description		RETAIL
Empulse M90 Option			\$
o 139MW94	Empulse M90 Wheels M90 Wheels, Dual Controller, T	wo M90 Batteries, Mounting Hardware, Quick Release Axles, and Device Charger	\$ 7,950.00
Options and Accessories			\$
o 139AG44	Additional Batteries	Includes two interchangeable M90 batteries, usable with either wheel.	\$ 500.00
o 139BJ39	Docking Charger	Includes docking charger, capable of charging two M90 batteries.	\$ 450.00
o 139JM50	Swing Away Mount Left Hand	Left hand swing away mount for M90 Dual Contoller.	\$ 175.00
o 139JM51	Swing Away Mount Right Hand	Right hand swing away mount for M90 Dual Contoller.	\$ 175.00
o 139NP02	Wrist support	Additional wrist and hand support for the M90 Dual Controller.	\$ 176.00
o 139DI56	Thigh Strap	Includes thigh strap and hardware, to mount the Dual Controller to the riders leg.	\$ 40.00
o 1399A33	Docking Tray	Docking tray for the M90 Dual Controller	\$ 30.00



WH	WHEELCHAIR										
а	OPTION	DESCRIPTION DESCRIPTION									
Sun	rise Medical Wheelchair Mo	dels - Specific mounting hardware	will be pi	rovide for the	chair selected						\$
Select the specific model of the wheelchair the M90 Wheels will be installed on.											
Quickie Rigid Chairs NC											
0	1399V58	Quickie Nitrum	О	1399V56	Quickie 5R		О	1399V59	Quickie Qri		
0	1399V80	Quickie GP		RGK -The E	RGK -The Empulse M90 works with select RGK models. Contact your local RGK rep		l RGK representativ	e for detai	ls.		
0	1399V66	Quickie 7R									
Qui	ckie Folding Chairs										NC
0	1399V69	Quickie QS5X	О	1399V68	Quickie 2		О	1399V60	Quickie Qxi		
0	1399V77	Quickie Qx	0	1399V78	Quickie 2 Lite	ickie 2 Lite					
Quickie Tilt In Space Chairs											NC
0	1399V05	Quickie Iris									
Zippie Pediatric - The Empulse M90 is only avalible in a 24" (540) wheel.											
0	4000104		0	1399V64 Zippie Zone		Zipi		Zippie X'Cape	ippie X'Cape Coming Soon		
0	1399V33	Zippie Iris									
Oth	Other Wheelchair Models - Specific mounting hardware will be provide for the chair selected \$										
		ir the M90 Wheels by Empluse will									
Sun	Sunrise Medical cannot guarantee fitment on all configurations of the below chairs, please contact your dealer for guidance.										
TiLi	ite Chairs - Coming Soon										NC
	TiLite ZR	TiLite TRA	TiL	Lite X							
	TiLite ZRA TiLite Aero Z		TiL	ite Aero T							
TiLite TR TiLite 2GX		TiLite 2GX	TiL	ite AeroX							
Ki Mobility Chairs - Coming soon										NC	
	Ki Mobility Catalyst ST	Ki Mobility Liberty FT									
	Ki Mobility Catalyst 5	Ki Mobility Focus CR									
Mot	tion Composites - Coming	soon									NC
	Apex A Helio C2		He	lio A6							
	Veloce	Helio A7									
Wh	eelchair Wheel Size										\$
) Wheels are only offered in										
M90) Wheels are only compatible	le with 0° of camber									
0	139WS6			2	4" Wheel Size (540 W	heels)					NC
Har	ndrim		\$								\$
Sele	ect only one Handrim option	(24" Wheels Only)									
О	139HR45	Powder Coated Aluminum Black	N	С		O 139HR39)	E	llipse 3R	\$	500.00
0	139HR4	Plastic Coated	\$ 1	129.00		O 139HR22		Natural F	it - Standard Grip	\$	549.00
0	139HR7	Omit Handrim	N	С							
Tire			\$	5							\$
Sel	ect only one Tire option (24"	Wheels Only)									
0	139RT21	Pneumatic Schwalbe Marathon Plus	N	С		o 139RT25		Prim	o Orion Poly		NC